

U-[¹³C₁₅]-DEOXYNIVALENOL IN ACETONITRILE

1. General information

This document is designed and the certified value(s) and uncertainty(ies) are determined in accordance with ISO Guide 31 [1] and Eurachem / CITAC Guides [2,3].

2. Description of the Reference Material (RM)

Name: U-[¹³C₁₅]-Deoxynivalenol in Acetonitrile
 CAS number: 51481-10-8
 Catalog number: B-MYC0330
 Lot #: I16194D
 Certificate version: 1
 Expiry date: 10.11.2017
 Starting material: U-[¹³C₁₅]-Deoxynivalenol, Lot #IS14472D, Romer Labs Diagnostic GmbH
 Physical description of RM: Solution of U-[¹³C₁₅]-Deoxynivalenol in Acetonitrile
 Packaging and amount of RM: Amber Certan® glass ampoules, solution of 1.2 mL
 Name and address of the manufacturer: Romer Labs Diagnostic GmbH
 Technopark 1, 3430 Tulln, Austria
 www.romerlabs.com
 Name and address of the supplier: LGC Standards GmbH
 Mercatorstraße 51, 46485 Wesel, Germany
 Tel +49(0)2 81 98 87 0, Fax +49(0)2 81/98 87 199
 www.lgcstandards.com

2.1 Intended use of the RM

- for laboratory use only
- internal standard

2.2 Instruction for the correct use of the RM

The ampoules should be stored at 2-8°C in a dark place. Before usage of the RM, the ampoules should be allowed to warm to room temperature. The recommended minimum sub-sample amount for all kinds of application is 100 µL. The expiry date of this RM is based on the current knowledge and holds only for proper storage conditions in the originally closed flasks/packages.

2.3 Hazardous situation

The normal laboratory safety precautions should be observed when working with this RM. Further details for the handling of this RM are available as safety data sheet.

Hazardous Ingredients	Concentration in %	Pictograms	Signal word	Hazard statement(s)
Acetonitrile	> 99.9		Danger	H225, H302, H312, H319, H332

3. Certified values and their uncertainties

U-[¹³ C ₁₅]-Deoxynivalenol in acetonitrile	
Compound	Mass concentration ^a

	Certified value ^b	Uncertainty ^c
U-[¹³ C ₁₅]-Deoxynivalenol, 99.0 atom % ¹³ C	25.4 µg/mL	± 0.2 µg/mL

^a Values are based on preparation data and confirmed experimentally by HPLC-UV
^b Mass concentration based on weighed amount, purity and dilution step
^c Expanded uncertainty U (k = 2) of the value u_c according to GUM [4]

3.1 Calculation of uncertainty

The uncertainty of the calibrant solution was calculated on the basis of preparation [5].

Uncertainty components	Description	Standard uncertainty (u)	
Purity (P) of solid U-[¹³ C ₁₅]-Deoxynivalenol, 99.0 atom % ¹³ C	P = 99.5 ± 0.5 %	u (P) = 0.3 %	a
Weighing procedure weighed sample: m _{ws} = 12.764 mg	U(m) = 0.0000008g + 1.21 * 10 ⁻⁵ * m _{Toxin} u(m) = U(m)/2	u (m) = 0.0005 mg	b
Dilution procedure volumetric flask: V _f = 500 mL	calibration: 500 mL ± 0.25 mL repeatability: 0.1 mL volume expansion solvent	u (cal) = 0.1 mL u (rep) = 0.1 mL u (Vol. exp.) = 1.2 mL u (V) = 1.2 mL	c d e f

^a Maximum tolerance of purity (rectangular distribution) was divided by $\sqrt{3}$

^b Calculation of this u-value is based upon the uncertainty formula for the weighed amount as given in the calibration report from annual balance calibration

^c A triangular distribution (division by $\sqrt{6}$) was chosen for the calculation of u (cal)

^d Based on a series of ten fill and weigh experiments on a typical 500 mL flask; the value was used directly as a standard deviation

^e Based on the density of 0.7857 g/cm³ at temperature T = 20°C and a maximum temperature variation of ± 3°C, of volume expansion, relative volume expansion coefficient of acetonitrile is 1370 * 10⁻⁶/°C [6], volume expansion term (rectangular distribution) was divided by $\sqrt{3}$

^f The three contributions are combined to give the u (V) = $\sqrt{u(\text{cal})^2 + u(\text{rep})^2 + u(\text{Vol. exp.})^2}$

Calculation of the combined uncertainty u_c and the expanded standard uncertainty U

$$c_{\text{Toxin}} = \frac{10' m_{\text{ws}}' P}{V_f} = \frac{10' 12.764' 99.5}{500} = 25.4 \text{ mg / L}$$

$$\frac{u_c(c_{\text{Toxin}})}{c_{\text{Toxin}}} = \sqrt{\frac{\dot{e}u(P)\dot{u}}{\hat{e}P\dot{u}} + \frac{\dot{e}u(m)\dot{u}}{\hat{e}m_{\text{ws}}\dot{u}} + \frac{\dot{e}u(V)\dot{u}}{\hat{e}V_f\dot{u}}} = \sqrt{\frac{\dot{e}0.3\dot{u}}{\hat{e}99.5\dot{u}} + \frac{\dot{e}0.0005\dot{u}}{\hat{e}12.764\dot{u}} + \frac{\dot{e}1.2\dot{u}}{\hat{e}500\dot{u}}} = 0.004$$


$$u_c(c_{\text{Toxin}}) = c_{\text{Toxin}}' 0.004 = 25.4' 0.004 = 0.1 \text{ mg / L}$$

Calculation of expanded standard uncertainty U using a coverage factor k = 2

$$U(c_{\text{Toxin}}) = u_c(c_{\text{Toxin}})' 2 = 0.1' 2 \gg 0.2 \text{ mg / L} = 0.2 \text{ mg / mL}$$

4. Isotopic enrichment and isotope pattern

Isotope pattern ^a	
Compound	Isotopic distribution
[¹³ C ₁₅]-Deoxynivalenol	86.5 %



[¹³ C ₁₄]-Deoxynivalenol	11.5 %	Figure 1: Enhanced resolution scan (250 amu/s) of U-[¹³ C ₁₅]-Deoxynivalenol for determination of isotope pattern
[¹³ C ₁₃]-Deoxynivalenol	2.1 %	
Calculated isotopic enrichment level ^a : 99.0 atom % ¹³C		
^a Approximation based on LC-MS/MS data		

5. Discussion of traceability

This calibrant is certified on the basis of gravimetric preparation [5]. Thus the certified value (mass concentration of U-[¹³C₁₅]-Deoxynivalenol, 99.0 atom % ¹³C) is based on the weighed amount of the starting material and is therefore traceable to the stated purity of the solid raw material. High purity material represents a practical realization of concentration units, through conversion of mass to molar quantity.

6. Confirmation of certified value by HPLC-UV

The certified concentration of U-[¹³C₁₅]-Deoxynivalenol, 99.0 atom % ¹³C) of the gravimetric prepared solution was confirmed by HPLC-UV against an independently prepared reference batch of unlabelled Deoxynivalenol calibrant Lot# L16123A.

column	Phenomenex Luna C18(2), 250x3.00 mm, 5μ	
injection volume	100 μL sample	
solvent A	acetonitrile / water (10 / 90)	
solvent B	acetonitrile	
flow rate	0.5 mL / min	
gradient	time in minutes (min)	% solvent B
	0 – 3	0
	3 – 6	0 - 11
	6 – 20	11
	20 – 21	11 - 0
	21 – 30	0
DAD settings	218 nm	
sample dilution	1:5 with acetonitrile / water (10 / 90)	




Figure 1: HPLC-UV chromatogram of U-[¹³C₁₅]-Deoxynivalenol calibrant Lot #16194D

	time [min]	area	concentration ^a [μg/mL]
U-[¹³ C ₁₅]-Deoxynivalenol	9.373	19.945	25.3 ± 0.8

^a Mean of 6 replicate measurements against reference batch, confidence interval with P = 95 %

7. Further information

The purchaser must determine the suitability of this product for its particular use. LGC Standards GmbH makes no warranty of any kind, express or implied, other than its products meet all quality control standards set by LGC Standards GmbH. We do not guarantee that the product can be used for a special application.

approved for release by: Dipl. Chem. Ulrich Berger

date: 13.05.2016

This document has been computer generated and is valid without a signature.

References:

- [1] ISO Guide **31**, 1-7, (2000), *"Reference Materials – Contents of Certificates and Labels"*
- [2] Eurachem / CITAC Guide, 1-37, (2003), *"Traceability in Chemical Measurement"*
- [3] Eurachem / CITAC Guide CG-4, 1-133, (QUAM:2012.P1), *"Quantifying Uncertainty in Analytical Measurement", 3rd Ed*
- [4] International Organization for Standardization (ISO), (2008), *"Guide to the expression of uncertainty in measurement"*, (GUM 1995 with minor corrections) 1st Ed. Geneva, Switzerland
- [5] R.D. Josephs, R. Krska, S. MacDonald, P. Wilson, H. Pettersson, J. AOAC Int. **86**, 50-60, (2003), *"Preparation of a Calibrant as Certified Reference Material for Determination of the Fusarium Mycotoxin Zearalenone"*
- [6] E.W. Flick, (1998), *"Industrial Solvents Handbook"*, 5th Ed., Noyes Data Corp. Westwood NJ